

mpgun®

m e s o p e n g u n



S1500



mpgun[®] S1500
mesopengun




Pen&Gun Performance

- 1500 stroke/min in Pen Mode
- 700 stroke/min Gun Mode
- 0.5 - 10 mm adjustable depth
- High penetration speed
- Adjustable frequency, pressure and depth parameters




User Friendly Menu

- 10 different recordable protocol memory



Safety

- Depth adjustment with 0.1 mm precision
- Resets at each start-up
- CE Class II a certified gun and needle kit



Special Design

- Hair and body application apparatuses (disposable)
- Special needle kit (disposable)
27Gx13mm, 30Gx13mm
- Easy to carry



Patient&Practitioner Comfort

- Painless application
- Saving on time and product





istanböl



BIOPARK[®]
MEDICAL

BIOPARK MEDİKAL SAN.TİC.A.Ş.
Sima İş Merkezi Fatih Cd. No. 18 K:3 D:10 34303
Halkalı / Küçükçekmece /İstanbul
T : +90 212 548 50 94 - F : +90 212 548 50 94
WhatsApp: +90 538 271 80 70
E : info@bioparkmedical.com